SPECIAL FEATURES & PRODUCT SPECIFICATIONS

Hydrophobic Filter Membrane
To facilitate blood flashback
and prevent any spillage of
blood during cannulation and
blood flashback.

Injection Port
Injection port enables
convenient injection of
medicine or during
infusion.

<u>Cannula</u> Ultra sharp cannula allows smooth needle

insertion.

Catheter Tip
Well formed catheter tip
provides an ultra-smooth
penetration, reduced
catheter peel back
and less pain and
trauma to patient.

Catheter Tube
Flexible and biocompatible
material ensures long hours
comfort and safety. Thin wall
allows maximum flow rate.

Flash Back Chamber
Allows monitoring of
instantaneous flash of
blood after venipuncture
confirming correct
catheter placement.

<u>Wing</u>
Flexible and foldable wings
for easy gripping and safe
fixing

14G	2.0	270	45	FEP	50/500
16G	1.7	180	45	FEP	50/500
18G	1.3	80	45	FEP	50/500
20G	1.1	54	32	FEP	50/500
22G	0.9	33	25	FEP	50/500
24G	0.7	20	19	FEP	50/500
26G	0.6	15	19	FEP	50/500



JMS SINGAPORE PTE LTD

440 Ang Mo Kio Industrial Park 1, Singapore 569620 TEL: (65) 64571144 FAX: (65) 64599564

www.jmss.com.sg

DISTRIBUTOR